
	<b>Document Information</b>	
	<b>File Name:</b> Raffaello AWS IHE Integration Statement v1.00.doc	<b>Version:</b> 1.00
	<b>Owner:</b> IMS s.r.l.	<b>Date:</b> 07-Feb-06

# Raffaello AWS IHE Integration Statement




Copyright © 2003-2006 by I.M.S. s.r.l.

	<b>Document Information</b>	
	<b>File Name:</b> Raffaello AWS IHE Integration Statement v1.00.doc	<b>Version:</b> 1.00
	<b>Owner:</b> IMS s.r.l.	<b>Date:</b> 07-Feb-06


## DOCUMENT VERSIONS

Version	Date	Author	Changes
1.00	07-Feb-06	IMS s.r.l.	First Version (Raffaello 1.1.3.7)

	<b>Document Information</b>	
	<b>File Name:</b> Raffaello AWS IHE Integration Statement v1.00.doc	<b>Version:</b> 1.00
	<b>Owner:</b> IMS s.r.l.	<b>Date:</b> 07-Feb-06

## TABLE OF CONTENTS

<b>0</b>	<b>INTRODUCTION .....</b>	<b>4</b>
<b>1</b>	<b>IHE INTEGRATION STATEMENT.....</b>	<b>5</b>


	<b>Document Information</b>	
	<b>File Name:</b> Raffaello AWS IHE Integration Statement v1.00.doc	<b>Version:</b> 1.00
	<b>Owner:</b> IMS s.r.l.	<b>Date:</b> 07-Feb-06

## 0 INTRODUCTION

The RAFFAELLO application is part of the following full field digital mammography systems manufactured from IMS:

- Giotto Image MD
- Giotto Image SD
- Giotto Image SDL

The RAFFAELLO software is used for both the Acquisition WorkStation (AWS) and the Review Workstation (RWS); this document contains only the IHE Integration Statement when the Raffaello software is configured for the AWS.

	<b>Document Information</b>	
	<b>File Name:</b> Raffaello AWS IHE Integration Statement v1.00.doc	<b>Version:</b> 1.00
	<b>Owner:</b> IMS s.r.l.	<b>Date:</b> 07-Feb-06

# 1 IHE INTEGRATION STATEMENT

<b>IHE Integration Statement</b>		<b>Date</b>	07 Feb 2006
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
IMS s.r.l.	Raffaello (configured for AWS)	1.1.3.7	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Consistent Presentation of Images	Print Composer	None	
Portable Data for Imaging	Portable Media Creator	None	
	Portable Media Importer	None	
	Image Display	None	
	Print Composer	None	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.imsitaly.com/IHE">http://www.imsitaly.com/IHE</a>			
<b>Links to Standards Conformance Statements for the implementation</b>			
<b>DICOM</b>	<a href="http://www.imsitaly.com/DICOM">http://www.imsitaly.com/DICOM</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	