
	Document Information	
	File Name: RaffaelloAWS_IHE_IntegrationStatement_v1.03.odt	Version: 1.03
	Owner: IMS s.r.l.	Date: 03-Nov-09

Raffaello AWS IHE Integration Statement

Copyright© 2003-2009 by I.M.S. s.r.l.

	Document Information	
	File Name: RaffaelloAWS_IHE_IntegrationStatement_v1.03.odt	Version: 1.03
	Owner: IMS s.r.l.	Date: 03-Nov-09

DOCUMENT VERSIONS

Version	Date	Author	Changes
1.00	07-Feb-06	IMS s.r.l.	First Version (Raffaello 1.1.3.7)
1.01	21-Feb-07	IMS s.r.l.	Updated for Raffaello 1.2.0.x
1.02	16-Oct-08	IMS s.r.l.	Updated for Raffaello 1.2.4.x
1.03	03-Nov-09	IMS s.r.l.	Updated for Raffaello 1.3.1.x



	Document Information	
	File Name: RaffaelloAWS_IHE_IntegrationStatement_v1.03.odt	Version: 1.03
	Owner: IMS s.r.l.	Date: 03-Nov-09

TABLE OF CONTENTS

0 INTRODUCTION.....4

1 IHE INTEGRATION STATEMENT.....5

	Document Information	
	File Name: RaffaelloAWS_IHE_IntegrationStatement_v1.03.odt	Version: 1.03
	Owner: IMS s.r.l.	Date: 03-Nov-09


0 INTRODUCTION

The RAFFAELLO application is part of the following full field digital mammography systems manufactured from IMS:

- **Giotto Image MD**
- **Giotto Image SD**
- **Giotto Image SDL**
- **Giotto Image 3D**
- **Giotto Image 3DL**

The RAFFAELLO application is also installed with the biopsy option for the analog mammography system Giotto Image (**Biopsy Digit-AM**) and with the **ConfirmX** system (radiologic equipment to image mammographic specimen).

The RAFFAELLO software is used for both the Acquisition WorkStation (AWS) and the Review Workstation (RWS); this document contains the IHE Integration Statement when the Raffaello software is configured for the AWS (this document can be used also in case Raffaello is installed as part of the Biopsy Digit-AM or on the ConfirmX).

	Document Information	
	File Name: RaffaelloAWS_IHE_IntegrationStatement_v1.03.odt	Version: 1.03
	Owner: IMS s.r.l.	Date: 03-Nov-09

1 IHE INTEGRATION STATEMENT

IHE Integration Statement		Date	03 Nov 2009
Vendor	Product Name	Version	
IMS s.r.l.	Raffaello (configured for AWS)	1.3.1.x	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
	Image Display	None	
Patient Information Reconciliation	Acquisition Modality	None	
Mammography Image	Acquisition Modality	Partial View	
	Print Composer	None	
Evidence Documents	Image Display	None	
Consistent Time	Time Client	None	
Consistent Presentation of Images	Print Composer	None	
Portable Data for Imaging	Portable Media Creator	None	
	Portable Media Importer	None	
	Image Display	None	
	Print Composer	None	
Internet address for vendor's IHE information: http://www.imsitaly.com/IHE/			
Links to Standards Conformance Statements for the implementation			
DICOM	http://www.imsitaly.com/DICOM/		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: http://www.ihe-j.org	